

Eurosafe 170

The XL automatic cleaning



Minimalist and compact, the Eurosafe 170 thermosinfectant meets the needs of large clinics. The large 170-litre chamber can hold up to 8 DIN 1/1 cassettes, equal to 4 loads of a 24-litre autoclave. The glass door allows the instruments to be monitored during the cycle and the pull-down opening serves as a convenient loading surface. It allows up to 20 programmes to be set and 3 favourite cycles to be combined with the buttons on the sides of the display, so that the desired wash can be selected immediately. The drying system is switched on and off using a convenient button. Eurosafe 170 monitors the correct regeneration of the integrated softener resins and warns the user of salt shortages via a pop-up message on the display. The 3.5" display shows the real-time temperature and cycle time. Efficiency and longevity are guaranteed by the injection washing and drying system with independent sprayers for each level and by a water filtering system that preserves the performance of the pump over time. In the event of malfunctions and leaks, the aquastop interrupts the water flow. The automatic dispenser independently calibrates and injects detergent liquids independently and without waste. Eurosafe 170 guarantees complete traceability: it validates and keeps track of all cycles performed and allows data to be saved on a USB stick, transferred to a PC via a data logger or printed out via an external printer (both optional). Made in Italy with durable AISI 316L stainless steel.



Classification

Medical Device Class IIb



Dimensions

60 x 84 x 63 cm (wxhxd)



Weight

90 Kg



Supply voltage

Single-phase 220 V / 50Hz
Three-phase 380 V / 50Hz

Power consumption

5600 W



Connectivity

Includes: USB and
E-Memory software
Optional: Datalogger



CE Ref

DS 500 CDL



Spacious

The 170-litre chamber can accommodate up to 4 full loads of a 24-litre autoclave. The numerous accessories available allow the best handling of handpieces, loose instruments and up to 8 DIN 1/1 cassettes.



Effective

The forced and filtered hot air system ensures complete internal and external drying. Eurosafe 170 optimises the cycle time, and energy and detergent consumption.



Powerful

The double wash pump with separate circuits ensures high water flow on both levels, while the two independent sprayers guarantee capillary action.



Full options

Eurosafe 170 facilitates load management and allows you to optimise thermal-disinfection cycles thanks to the wide range of accessories included: the upper and lower trolley, the insert for vertically-positioned instruments, the DIN 1/1 basket, and the 2 8-slot tray inserts.

