




检验报告



Information Provided by Client	Applicant	Changzhou Jinliyan Medical Devices Co., Ltd. Contact --- No.6, West Xihe Road, Sanhekou, Zhenglu Town, Tianning District, 213115 Changzhou City, P. R. China Tel. ---	
	Manufacturer	Changzhou Jinliyan Medical Devices Co., Ltd. No.6, West Xihe Road, Sanhekou, Zhenglu Town, Tianning District, 213115 Changzhou City, P. R. China	
	Information of Submitted Sample	Sample Name	Disposable Medical Face Mask (LOT: 20200508) Trademark -----
		Sample Count	62Piece(s)
	Size	Ear Loop 175mm×95mm Colour -----	
	Quality Grade	----- Safety Category -----	
	Style No. or Order No.	-----	
Test Part Description	1# Blue Mask		
Test Type	Commission Test	Date of Submission	2020-05-16 Date of Checking 2020-06-01
Test Standards	BS EN 14683:2019 Medical face masks-Requirements and test methods		
Conclusion	Test results and compliance refer to next page(s). Stamp of Inspection Unit		
Non-standard Check Methods	-----		
Uncertainty of Result	-----		
Remarks	----- Face side up		



Approver

方倩 路妍

Checker

张涛

Editor

李刚



检验报告

CTTT-WT20011607

Page 2 of 2

Test Items	Description	Unit	Standard Requirement	Results	Conclusions	Test Method/Remarks
1#						
Splash Resistance Pressure (Type IIR)	---	kPa	≥16.0	≥16.0	Pass	ISO 22609:2004
Differential Pressure (Type IIR)	---	Pa/cm ²	<60	41	Pass	BS EN 14683:2019
Bacterial Filtration Efficiency (BFE) (Type IIR)	---	%	≥98	1: >99.9 2: >99.9 3: 99.7 4: 99.8 5: 99.8	Pass	BS EN 14683:2019
Microbial Cleanliness* (Type IIR)	---	CFU/g	≤30	14	Pass	BS EN 14683:2019 EN ISO 11737-1:2018

* Indicates the item is unauthorized by CMA

Blank Part Below